

Metrohm *tiBase* – intelligent management of titration data

***tiBase* brings database functionality to Metrohm's Titrino plus and Titrando titrators. The new software enables the straightforward transfer of titration determinations via USB memory stick to a PC where the data can be stored, administrated, reprocessed and protected against any manipulation.**



Database functionality, previously only available in Metrohm's high-end range of instruments, is now also available in Metrohm's entry level range, namely the Titrino plus titrators. ***tiBase*** comes with a multitude of data management functions such as

- archiving
- easy data import and export
- search and filter functions
- reprocessing
- freely definable report formats

tiBase provides for extensive data security. Not only determinations, but all settings, methods and user information are stored in the database and protected against wilful manipulation or loss. ***tiBase*** Multi is the client-server version of the new software and enables central data administration by way of client-server option.